


DB Plus COMPOSITE POSTERIOR TEETH
DB Plus COMPOSITE POSTÉRIEURES

20° COMPOSITE POSTERIOR TEETH
20° COMPOSITE POSTÉRIEURES



20° COMPOSITE POSTERIOR TEETH
20° COMPOSITE POSTÉRIEURES



30° COMPOSITE POSTERIOR TEETH
30° COMPOSITE POSTÉRIEURES



30° COMPOSITE POSTERIOR TEETH
30° COMPOSITE POSTÉRIEURES



db plus COMBINATION CHART

MAXILLARY ANTERIOR MOULD	ARTICULATES WITH MANDIBULAR MOULD	COMBINES WITH 0° POSTERIOR MOULD	COMBINES WITH F-10 POSTERIOR MOULD	COMBINES WITH 20° POSTERIOR MOULD	COMBINES WITH 30° POSTERIOR MOULD
O40	L1	29	F29	29M or 29L	S30/L30/31
O45	L1	29	F29	29M or 29L	S30/L30/31
O52	L1	29	F29	29M or 29L	S30/L30/31
O55	L8	29	F29	29M or 29L	S30/L30/31
O65	L9	31	F31	31M or 31L	S32/L32/33
O75	L13	31	F31	31M or 31L	S32/L32/33
O80	L13	33	F33	33M	33/X34
O89	L13	33	F33	33M	33/X34
O99	L16	33	F33	33M	33/X34
T40	L1	29	F29	29M or 29L	S30/L30/31
T43	L3	29	F29	29M or 29L	S30/L30/31
T53	L3	29	F29	29M or 29L	S30/L30/31
T54	L5	29	F29	29M or 29L	S30/L30/31
T63	L6	31	F31	31M or 31L	S32/L32/33
T76	L11	31	F31	31M or 31L	S32/L32/33
T81	L9	33	F33	33M	33/X34
T90	L15	33	F33	33M	33/X34
R40	L1	29	F29	29M or 29L	S30/L30/31
R41	L1	29	F29	29M or 29L	S30/L30/31
R58	L3	31	F31	31M or 31L	S32/L32/33
R61	L5	31	F31	31M or 31L	S32/L32/33
R69	L9	33	F33	33M	33/X34
R78	L15	33	F33	33M	33/X34
R88	L13	33	F33	33M	33/X34
R91	L15	33	F33	33M	33/X34
S45	L1	29	F29	29M or 29L	S30/L30/31
S56	L5	29	F29	29M or 29L	S30/L30/31
S65	L7	31	F31	31M or 31L	S32/L32/33
S67	L7	31	F31	31M or 31L	S32/L32/33
S71	L11	31	F31	31M or 31L	S32/L32/33
S79	L13	31	F31	31M or 31L	S32/L32/33
S83	L13	33	F33	33M	33/X34
S98	L13	33	F33	33M	33/X34

myerson.

MANUFACTURER
Myerson Company Limited
3 trinity avenue
Laventille, Trinidad & Tobago
phone 868 623 1007


DISTRIBUTED BY
Myerson LLC
5106 North Ravenswood
Chicago IL 60640-2713
USA
phone 800 423 2683

EC REP
Obelis s.a
Bd. Général Wahis 53
1030 Brussel, Belgium
phone +(32) 2.732.59.54
fax +(32) 2.732.60.03
email: mail@obelis.net

Preformed plastic denture teeth for use in making temporary or removable dentures.
Documented as per ISO 22112:2017
2020-01-21

Make great removables
with Myerson teeth
& Myerson's advanced
thermoplastic materials

DURAFLEX[™]
PARTIALS



VISICLEAR[®]
CLASPS & FRAMEWORKS



DURACETAL[®]
COSMETIC APPLIANCES



ema[®]
from
myerson.

With 4 strap strengths and
9 strap lengths, a custom
EMA treats snoring and
obstructive sleep apnea
safely and effectively, giving
you the sleep you need to
lead a full, healthy life.



ema[®] is a trademark of frantz design inc.
US Patent Numbers: 6,109,265; D611,153 & Patent Pending

myerson DB Plus
mould chart

myerson[®]

We make the teeth
that make the smile.

DB Plus RESIN ANTERIOR TEETH • MAXILLARY
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE SUPÉRIEUR

OVOID EN FORME DE PELLE			TAPERING TRIANGULAIRE		
O40		39.7 8.5 7.0 (L1)	T40		40.3 9.5 7.3 (L1)
O45		41.5 9.4 7.9 (L1)	T43		41.4 9.6 7.6 (L3)
O52		42.4 8.7 7.8 (L1)	T53		42.8 9.9 7.9 (L3)
O55		43.1 9.4 8.0 (L8)	T54		43.1 9.9 7.7 (L5)
O65		45.0 9.9 8.3 (L9)	T63		44.5 10.5 8.0 (L6)
O75		46.8 11.1 8.7 (L13)	T76		46.9 10.3 8.5 (L11)
O80		47.3 10.9 8.8 (L13)	T81		48.0 11.1 8.8 (L9)
O89		49.0 11.3 9.6 (L13)	T90		49.4 12.4 9.3 (L15)
O99		51.5 12.1 9.4 (L16)			

TOOTH SHAPE KEY		LÉGENDE DE FORMES DES DENTS	
O	T	S	R
OVOID EN FORME DE PELLE	TAPERING TRIANGULAIRE	SQUARE CARRÉ	RECTANGULAR RECTANGULAIRE



DB Plus RESIN ANTERIOR TEETH • MAXILLARY
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE SUPÉRIEUR

RECTANGULAR RECTANGULAIRE			SQUARE CARRÉ		
R40		39.5 9.6 7.0 (L1)	S45		41.5 8.7 7.4 (L1)
R41		40.5 9.7 7.2 (L1)	S56		43.2 9.2 7.9 (L5)
R58		43.5 9.9 7.9 (L3)	S65		45.0 10.6 8.1 (L7)
R61		43.8 10.0 8.1 (L5)	S67		45.4 9.6 8.2 (L7)
R69		45.8 10.4 8.3 (L9)	S71		46.2 10.1 8.4 (L11)
R78		47.0 11.4 8.6 (L15)	S79		47.0 9.2 9.1 (L13)
R88		49.0 10.4 9.1 (L13)	S83		48.2 10.2 8.8 (L13)
R91		50.3 12.9 9.3 (L15)	S98		51.4 10.5 9.5 (L13)

DB Plus RESIN ANTERIOR TEETH • MANDIBULAR

L1		31.0 8.5 4.4
L3		32.0 9.1 4.6

DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE INFÉRIEUR

L5		32.9 8.5 4.6
L6		33.1 9.8 5.0

DB Plus RESIN ANTERIOR TEETH • MANDIBULAR
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE INFÉRIEUR


L7		33.6 9.9 4.9
L8		33.7 8.5 4.6
L9		34.7 9.8 5.1
L11		35.0 9.0 5.2
L13		36.4 9.0 5.5
L15		36.7 11.1 5.6
L16		40.0 10.5 5.9

DB Plus COMPOSITE POSTERIOR TEETH
DB Plus COMPOSITE POSTÉRIEURES

0° COMPOSITE POSTERIOR TEETH • 0° COMPOSITE POSTÉRIEURES

29		29.0 30.5 7.10 7.85	
31		30.5 32.3 7.55 8.50	
33		32.8 34.8 8.02 8.74	

F10 COMPOSITE POSTERIOR TEETH • F10 COMPOSITE POSTÉRIEURES

F29		28.6 31.6 8.36 7.48	
F31		31.0 31.8 9.19 8.37	
F33		32.8 34.3 9.62 9.09	
F35		36.5 37.5 10.8 9.62	

myerson.