


KENSON F35

F35

Introducing F35, the new extra large mould size in Kenson F-10 Posterior Teeth. This natural looking, shallow cusp posterior is wider and longer than the popular F33 mould size, making F35 ideal for partial dentures. Available in 16 matching VITA shades and whiter white shade 51.




KENSON XL & X36 TEETH

XL Teeth

Upper central incisors only...over two millimeters wider and longer*. Available in all Kenson shades.

The new Kenson XL mould saves you the trouble and expense of making teeth from scratch or reaching for more expensive brands when you have a large anterior space to fill in an upper complete or partial denture



X36 Teeth


Upper bicuspids only... Kenson mould X36 is an ideal space filler for partial dentures; can also be used as a molar.



KENSON RESIN VITA SHADES


Kenson Teeth, the professional's choice for 50 years, are now available in 16 matching VITA shades that are True to Life.

Myerson continues to offer Kenson Teeth in the 10 original Kenson shades: 55, 61, 62, 65, 66, 67, 69, 77, 81, 87, plus New whiter white shade 51.



myerson Part Number T1SGV

KENSON POSTERIOR KEY



Kenson Resin COMBINATION CHART

MAXILLARY ANTERIOR MOULD	ARTICULATES WITH MANDIBULAR MOULD	COMBINES WITH 0° POSTERIOR MOULD	COMBINES WITH F-10 POSTERIOR MOULD	COMBINES WITH 20° POSTERIOR MOULD	COMBINES WITH 30° POSTERIOR MOULD
14	43	29	F29	29M or 29L	S30•L30•31
15	47	31	F31	31M or 31L	S32•L32•33
16	49	31	F31	31M or 31L	S32•L32•33
17	49	33	F33	33M	33•X34
27	39	29	F29	29M or 29L	S30•L30•31
28	46	31	F31	31M or 31L	S32•L32•33
29	48	31	F31	31M or 31L	S32•L32•33
P5	46	29	F29	29M or 29L	S30•L30•31
P6	48	31	F31	31M or 31L	S32•L32•33
P7	49	33	F33	33M	33•X34
S5	53	33	F33	33M	33•X34
S6	53	33	F33	33M	33•X34
T3	43	29	F29	29M or 29L	S30•L30•31
T4	46	29	F29	29M or 29L	S30•L30•31
T33	47	29	F29	29M or 29L	S30•L30•31
T34	47	31	F31	31M or 31L	S32•L32•33
T5	48	31	F31	31M or 31L	S32•L32•33
T6	48	33	F33	33M	L32•33

*0 and "BLOCK" Posterior teeth are interchangeable



Preformed plastic denture teeth for use in making temporary or removable dentures.

DISTRIBUTED BY
Myerson LLC
5106 North Ravenswood, Chicago, IL 60640-2713, United States
toll free phone 800 423 2683, phone 312 432 8200, fax 800 424 2928

myersontooth.com


EUROPEAN REPRESENTATIVE
Obelis s.a
Bd. Général Wahis 53
1030 Brussel, Belgium
phone +(32) 2.732.59.54
fax +(32) 2.732.60.03
email: mail@obelis.net

MANUFACTURED BY
Myerson Company Limited
3 Trinity Avenue
Laventille
Trinidad & Tobago
phone 868 623 1007

Documented as per ISO 22112:2017
Rev A 2020-08-27

Make great removables with Myerson teeth & Myerson's advanced thermoplastic materials

DURAFLEX PARTIALS



VISICLEAR CLEAR & TRANSPARENT



DURACETAL COSMETIC APPLIANCES



Kenson® Teeth
The Professional Choice Since 1955

ema
from
myerson.

With 4 strap strengths and 9 strap lengths, a custom EMA treats snoring and obstructive sleep apnea safely and effectively, giving you the sleep you need to lead a full, healthy life.



ema is a trademark of franz design inc.
US Patent Numbers: 6,109,200, D611,153 & Patent Pending

KENSON RESIN ANTERIOR TEETH • MAXILLARY


14 △		43.1 9.9 7.7 (43)	P7 □		48.2 10.2 8.8 (49)
15 □		45.0 10.6 8.1 (47)	S5 ○		51.5 12.1 9.4 (53)
16 ○		46.8 11.0 8.8 (49)	S6 □		51.4 10.5 9.5 (53)
17 ○		49.0 11.3 8.8 (49)	T3 △		41.7 9.7 7.5 (43)
27 △		41.5 8.7 7.4 (39)	T4 △		44.6 9.8 8.4 (46)
28 □		45.4 9.6 8.2 (46)	T33 △		44.0 10.2 8.3 (47)
29 □		46.9 10.3 8.5 (48)	T34 △		45.6 10.8 8.3 (47)
P5 ○		43.1 9.4 8.0 (46)	T5 □		47.0 9.2 9.1 (48)
P6 □		46.2 10.1 8.4 (48)	T6 □		47.9 10.2 8.5 (48)

TOOTH TYPE



KENSON RESIN ANTERIOR TEETH • MANDIBULAR

39		31.0 8.5 4.4	49		36.4 9.0 5.5
43		32.6 9.1 4.9	50		37.0 11.3 5.6
46		33.7 8.5 5.0	52		37.0 12.3 5.5
47		33.6 9.9 4.9	53		40.0 10.5 5.9
48		35.6 8.8 5.4	14 △		43.1 9.9 7.7 (43)



KENSON RESIN • 0° POSTERIOR TEETH

29		29.0 30.5 7.10 7.85
31		30.5 32.3 7.55 8.50
33		32.8 34.9 8.02 8.74

KENSON RESIN • F10 POSTERIOR TEETH

F29		28.6 31.6 8.36 7.48
F31		31.0 31.8 9.19 8.37
F33		32.8 34.3 9.62 9.09
F35		36.5 37.5 10.8 9.62

KENSON RESIN • 30° POSTERIOR TEETH

S30		28.5 29.4 9.12 8.53	L32		31.2 31.1 9.15 9.15
L30		28.8 30.2 9.10 9.15	33		31.8 35.7 9.60 9.82
31		31.3 33.9 9.50 10.15	X34		33.0 35.8 9.94 9.74
S32		31.6 32.0 9.54 9.09	X36		15.7 *Width of 2

KENSON RESIN • BLOCK POSTERIOR TEETH

29		28.7 30.4 7.10 7.85
31		30.5 32.4 7.30 8.50
33		33.0 34.5 8.00 8.70

KENSON RESIN • 20° POSTERIOR TEETH

29M/L		29.0-31.6/29.0-31.4 8.10-8.85/8.85-8.70
31M/L		31.4-31.6/31.6-34.8 8.10-9.10/9.55-9.25
33M		33.1 35.8 8.35 9.40