Make great removables with Myerson teeth & Myerson's advanced thermoplastic materials





DURA CETAL



myerson.

Preformed plastic denture teeth for use in making temporary or removable dentures.

DISTRIBUTED BY Myerson LLC 5106 North Ravenswood, Chicago, IL 60640-2713, United States toll free phone 800 423 2683, phone 312 432 8200, fax 800 424 2928

myersontooth.com



MANUFACTURED BY **Myerson Company Limited** 3 Trinity Avenue Laventille Trinidad & Tobago phone 868 623 1007

Documented as per ISO 22112:2017 Revision 2020-01-21



With **4** strap strengths and **9** strap lengths, a custom EMA treats snoring and obstructive sleep apnea safely and effectively, giving you the sleep you need to lead a full, healthy life.



Myerson Special Mould Chart

MARSON.

We make the teeth that make the smile.

myerson teeth mirror the beauty & strength of natural teeth.



Approximal stains Strong surface relief Enamel striation Darker neck of tooth





Myerson provides a premium denture solution with Myerson Special® teeth that meet the profession's highest standard in denture aesthetics. Myerson Special teeth come in an extensive selection of personalized moulds and all 16 A-D shades.


Myerson Special Aesthetic Advantages

- Natural fluorescence: consistent shade values under any light conditions.
- Natural translucency & surface luster: harmonious with adjacent natural teeth.
- Natural labial carvings: perfectly suited for partial dentures

Characterizing features are individually added by hand to simulate the appearance of natural dentition. Colors and special effects come from the tooth rather than on the surface, which means that the teeth can be modified without affecting the shade. The result is a look that will satisfy even the most discerning patient.



We make the teeth that make the smile.




stinguishable from natural teeth, than


Myerson Special anterior teeth combin Myerson DB Plus composite posterior in which our Durablend resin is reinfor with sub-micron silica, further enhanc

Composite filling Decalcification Translucent incisal areas Darker canine teeth.


COMBINATION CHART

		ARTICULATES	COMBINES WITH DB PLUS 0	COMBINES WITH	COMBINES WITH	COMBINES WITH
	MAXILLARY ANTERIOR	WITH Mandibular	POSTERIOR	DB PLUS F-10 POSTERIOR	DB PLUS 20 POSTERIOR	DB PLUS 30 POSTERIOR
	MOULD	MOULD	MOULD	MOULD	MOULD	MOULD
	A2	11	29	F29	29M or 29L	S30, L30 or 31
	A4	7	29	F29	29M or 29L	S30, L30 or 31
	A7	7	31	F31	31M or 31L	S32, L32 or 33
	A9	8	31	F31	31M or 31L	S32, L32 or 33
	A11	14	33	F33	33M	33 or X34
	A13	16	33	F33	33M	33 or X34
	A43	12	29	F29	29M or 29L	S30, L30 or 31
	A44	13	31	F31	31M or 31L	S32, L32 or 33
	A52	13	33	F33	33M	33 or X34
	A53	16	33	F33	33M	33 or X34
	A64	A71	33	F33	33M	33 or X34
	A72	A73	33	F33	33M	33 or X34
	YM3	11	29	F29	29M or 29L	S30, L30 or 31
·	YM5	8	31	F31	31M or 31L	S32, L32 or 33
	YW3	13	31	F31	31M or 31L	S32, L32 or 33
	YW5	14	31	F31	31M or 31L	S32, L32 or 33
·	PM2	11	29	F29	29M or 29L	S30, L30 or 31
	PM4	12	31	F31	31M or 31L	S32, L32 or 33
·	PM6	13	31	F31	31M or 31L	S32, L32 or 33
·	PM8	14	33	F33	33M	33 or X34
	PW2	13	31	F31	31M or 31L	S32, L32 or 33
	PW4	A73	33	F33	33M	33 or X34
	PW6	16	33	F33	33M	33 or X34
	PW8	16	33	F33	33M	33 or X34





C



10.8 9.62

20 DEGREE

